

b. Bacterial and Eukaryotic UTR Insertions

In general, RNA Mango aptamers can be easily inserted in the 5' and 3' UTR of a target mRNA:

i. For a 5'UTR insertion, use caution when inserting in this region. Ensure that the RNA Mango aptamer being inserted does not:

- interrupt any known existing structural RNA elements
- contain the start codon
- interrupt the Shine-Dalgarno sequence in bacteria
- disrupt eukaryotic ribosomal scanning

ii. For a 3'UTR insertion, add a Mango in any frame. Mango II aptamer arrays have been inserted immediately or shortly after the stop codon in a gene of interest to successfully track RNA:

- <https://www.nature.com/articles/s41467-020-14932-7>

c. ORF Insertions

For ORF insertions, Mango aptamers can be added in any frame and will result in the alteration of a local region of the protein sequence. You may use the alignments below to determine how these insertions can change the protein that will be translated. **Conveniently, the dye binding cores of Mango I, II, and III/IIIA10U aptamers do not contain stop codons in any frame.** To keep the aptamers and the remainder of the peptide sequence in frame, simply add nucleotides (+ 1 or + 2) as illustrated in the examples below.

	Reading Frame 1	Reading Frame 2	Reading Frame 3	(X dependent on stem chosen)							
Mango I	NNN NNN	CGA R	AGG R	GAC D	GGU G	GCG A	GAG E	AGG R	AGA R	GN'N' X	N'N'N' N' + 2
	NN NNN	NCG X	AAG K	GGA G	CGG R	UGC C	GGA G	GAG E	GAG E	AGN' X	N'N'N' N'N' + 2
	N NNN	NNC X	GAA E	GGG G	ACG T	GUG V	CGG R	AGA R	GGA G	GAG E	N'N'N' N'N'N' + 2
Mango II	NNN NNN	CGA R	AGG R	AGA R	GGA G	GAG E	GAA E	GAG E	GAG E	AGN' X	N'N'N' N'N' + 1
	NN NNN	NCG X	AAG K	GAG E	AGG R	AGA R	GGA G	AGA R	GGA G	GAG E	N'N'N' N'N'N' + 1
	N NNN	NNC X	GAA E	GGA G	GAG E	GAG E	AGG R	AAG K	AGG R	AGA R	GN'N' X
Mango III	NNN NNN	AGG R	AUU I	GGU G	AUG M	UGG W	UAU Y	AUN' X			N'N'N' N'N' + 1
	NN NNN	NAG X	GAU D	UGG W	UAU Y	GUG V	GUA V	UAU Y			N'N'N' N'N'N' + 1
	N NNN	NNA X	GGA G	UUG L	GUA V	UGU C	GGU G	AUA I	UN'N' X		N'N'N' N' + 1
Mango III (A10U)	NNN NNN	AGG R	UUU F	GGU G	AUG M	UGG W	UAU Y	AUN' X			N'N'N' N'N' + 1
	NN NNN	NAG X	GUU V	UGG W	UAU Y	GUG V	GUA V	UAU Y			N'N'N' N'N'N' + 1
	N NNN	NNA X	GGU G	UUG L	GUA V	UGU C	GGU G	AUA I	UN'N' X		N'N'N' N' + 1