

Service Report

Order 12345
Mycoplasma Detection
Date

Applied Biological Materials Inc.
1-3671 Viking Way,
Richmond, BC, CANADA
V6V 2J5

Service Description

Service Name:	Mycoplasma Detection
Service Summary:	PCR based mycoplasma detection. Mycoplasma contamination is indicated by the presence of bands at approximately 500bp.
Cat. Number:	C214
Sample(s) Information:	Sample 1
Deliverable:	Service Report

Test Results:

Test	Method	Results
Mycoplasma	PCR	Not Detected
Sterility	Direct Culture	Not Detected

PCR Amplification Gel:

Legend:

L: 100 bp DNA Ladder
NTC: Non-Template Control
Sample 1
P: Positive Control

