


GFP Control Lentivirus (Puro)



Cat. No.

LV095

Quantity

10 ml

Species	N/A
Accession Number	N/A
Vector	pLenti-III-CMV-GFP Control Vector
Vector Size	8934bp
Insert Size	N/A
Sequencing Primer	5'---CGC AAA TGG GCG GTA GGC GTG---3'
Tags	N/A
Selection Marker	Mammalian: Puromycin
Format	Lentivirus
Appearance	Liquid
Storage Buffer	Serum-free DMEM
Storage	1 year when stored at -80°C.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).