



Certificate of Analysis

Product Description

Product Name	Immortalized Hair Follicle Dermal Papilla Cells – SV40
Cat Number	T0500
Lot Number	HC0613
Source Organism	<i>Homo sapiens</i>
Source Organ	Hair Follicle
Viability	>70%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	June 27, 2013
QC Evaluation Medium	Prigrow III (TM003)

Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected

Cell Type Characterization

Gene	Cell Type	Specification		Results			
SFRP2	Hair Follicle Dermal Papilla	Positive	Pass	SFRP2	SMA- α	GAPDH	
Smooth Muscle Actin- α							
GAPDH	House-Keeping	Positive	Pass	Size	776bp	862bp	208bp

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

No. 8, 13520 Crestwood Place
 Richmond BC, Canada V6V 2G2
 T e l : 6 0 4 - 2 4 7 - 2 4 1 6
 F a x : 6 0 4 - 2 4 7 - 2 4 1 4
 w w w . a b m G o o d . c o m